Stehgeräte

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Instruction for use

Stand: 09/2022
(Rev. 4.0)

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6.5	Thoracic Pad	41
6.5.1	Setting the Thoracic Pad.....	41
6.5.2	Thoracic Lateral Pads.....	41
6.5.3	Adjusting the Width Setting of the Thoracic Lateral Pads.....	41
6.5.4	Fitting and adjusting the Head Support	41
6.6	Thorax-Support.....	42
6.7	General instructions for use when standing up	42
6.8	Operating instructions.....	44
6.9	Behavior in case of power failure.....	45
6.10	Operating the brake castors.....	45
7	Ambient conditions	46
7.1	Storage conditions.....	46
7.2	Operating conditions.....	46
8	Technical data.....	47
9	Used materials.....	48
10	Service and care	48
11	Service life of the product	48
12	Disinfection.....	48
12.1	Specifications of detergents and disinfectants	48
13	Operational faults and solutions	49
14	Recommended accessories	49
15	Maintenance.....	50
15.1	Legal basis.....	50
15.2	Maintenance intervals.....	50
15.3	Spare parts.....	50
15.4	Notes on documentation	51
16	Reuse	51
17	Disposal	51
17.1	Disposal of the device.....	51
17.2	Disposal of the electrical components	51
17.3	Disposal of the packaging.....	51
18	Declaration of Conformity.....	52

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List of figures

Figure 1: Exemplary type plate.....	34
Figure 2: Adjustment of the table height	39
Figure 3: Clamping lever for table height adjustment	39
Figure 4: Positioning the lateral pads.....	40
Figure 5: Adjustment of heel cups and toe fixation	40
Figure 6: Adjustment of the knee pads	41
Figure 7: Thoracic Pad	42
Figure 8: Thorax-support.....	42
Figure 9: Position of the patient support belt in the buttocks when the patient is in the wheelchair.	43
Figure 10: Hand control.....	43
Figure 11: Electric retractable belt	44
Figure 12: Exemplary castors in braked and unbraked condition.....	45

List of tables

Table 1: Used symbols.....	34
Table 2: Standards verification.....	35
Table 3: Storage conditions.....	46
Table 4: Operating conditions	46
Table 5: Technical data.....	47

Table 6: Operational faults and solutions	49
Table 7: Recommended accessories	49

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1 Foreword

Dear customer!



We would like to thank you for the trust you have placed in us and for purchasing our product. We have manufactured this medical product with great care.

Please read the instructions for use carefully before using the product for the first time and always keep them close at hand.

Not all conceivable uses of the device can be covered in these instructions for use. For further information or in the event of problems that are not described in sufficient detail in these instructions for use, please contact your specialist dealer or medical supply store.

2 General notes

2.1 Used symbols

	<p>This warning sign indicates all instructions that are important for safety. Non-observance can lead to accidents or injuries.</p>
	<p>Manufacturer - Indicates the manufacturer of the medical device according to EU Directives 2017/745. The symbol must appear in close proximity to the symbol, together with the name and address of the manufacturer (i.e. the person who places the medical device on the market)</p>
	<p>Conformity symbol according to 2017/745 of the Medical Devices Directive</p>
	<p>Medical Device - Shows the medical device provided by the manufacturer in accordance with EU Directives 2017/745</p>
	<p>Device type B according to IEC 601-1 (Special protection against electric shock)</p>
	<p>Device of protection class II, protective insulation</p>
	<p>Dispose of electrical components in accordance with the legal requirements. Do not dispose of in household waste!</p>
	<p>Date of manufacture - indicates the date when the medical device was manufactured.</p>
	<p>Part number - displays the manufacturer's part number so that the medical device can be identified.</p>

	<p>Serial number - displays the manufacturer's serial number so that a specific medical device can be identified.</p>
	<p>Distributor - indicates the company that distributes the medical device at the location.</p>
	<p>Temperature Limit - indicates the temperature limits to which the medical device can be safely exposed.</p>
	<p>Humidity, Limit - indicates the humidity range to which the medical device can be safely exposed.</p>
	<p>Air Pressure, Limit - indicates the range of air pressure to which the medical device can be safely exposed.</p>
	<p>Observe instruction for use or electronic instruction for use - indicates to the user that it is necessary to observe the instruction for use.</p>
	<p>Unique identifier of a medical device - displays a carrier containing information about a unique identifier of a medical device.</p>
	<p>Safe working load</p>
	<p>Max. patient weight</p>

	Minimum body dimensions/weights of the patient
--	--

Table 1: Used symbols

2.2 Type plate

The type plate is attached to the head of the trolley frame. The nameplate allows the product to be clearly identified.

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13.01.2022

Figure 1: Exemplary type plate

Figure 1 shows an exemplary type plate. For the exact specifications of your product, please refer to the attached type plate.

2.3 Standards verification

The following national and international norms (standards) are used in the design and verification of the product, labeling and instructions for use.

Standard	Title	Edition
DIN EN 60601-1-6	<i>Serviceability specification</i>	2010
EN 60601-1-2	<i>Electromagnetic compatibility</i>	2015
DIN EN ISO 10993	<i>Biological evaluation of medical devices - Part 1: Assessment and testing</i>	2010
DIN EN 1041	<i>Provision of information by the manufacturer of a medical device</i>	2008
DIN EN ISO 14971	<i>Medical devices - Application of risk management to medical devices</i>	2020

Table 2: Standards verification

3 Safety instructions

- You should read this instruction for use carefully (see Medical Devices Operator Ordinance under your national law) before using the standing frame. It contains important information for the safe and reliable use of the device. Keep the instruction for use for future reference.
- Safety, reliability and performance are guaranteed if the following instructions are observed and the device is used in an expert manner. As the operator, you must comply with the Medical Devices Operator Ordinance under your national law.
- The device may only be operated by medically trained and instructed personnel when standing training begins. Operators must be trained in the use of the device and their training and level of knowledge must allow an assessment of the risks for the patient. The trained personnel will take over the assessment for subsequent independent use by the user (patient), whereby an assistant must remain within call or sight.
- All clamp and screw connections must be checked regularly
- During use, 1 or 2 helpers must be available to assist the user
- Only use aids on a firm, level and horizontal surface.
- Do not use the device outdoors. The device may only be used in protected domestic and clinical areas.
- Leaning out far poses a risk of falling
- The medical device must not be used as a means of transport
- Park the standing frame only with locked wheels
- In the event of a malfunction, the standing frame must be immediately withdrawn from use
- Parts of the standing frame can catch fire if handled improperly
- when standing, the user must always wear sturdy shoes, barefoot use is not permitted. Unless expressly instructed by the therapist
- Ensure that children only have access to the standing frame under supervision.
- The standing frame should only be set up by authorized personnel.
- The installation-side fuse protection must not exceed 16A. Before connecting the device, please make sure that the voltage and frequency of your mains supply correspond to the specifications on the type plate.
- Ensure a level surface when selecting the location
- Provide a suitable floor covering if the standing frame has to be moved frequently. Carpets, rugs and loosely laid floor coverings can be damaged or make pushing difficult.
- Connect the mains plug firmly to the mains socket. Lay the power supply cable on the floor. Make sure that the standing frame (especially when moving it) does not rest on the cable with its rollers. The cable must not be routed through the mechanical parts of the base unit! (Danger of crushing)
- Damaged power cables can lead to life-threatening situations. These must be replaced immediately.
- Check the power cord for damage at regular intervals (weekly).
- Make sure that the electrical specifications of the device match the local conditions at the installation site.
- When the hand control are not in use, make sure that they are hanging on the standing frame to prevent incorrect operation that could cause damage.

4 General product description

4.1 Intended purpose

The standing device is used exclusively for standing training for people who are unable to walk or have difficulty walking for independent use with an assistant. Standing training is used in inpatient and outpatient rehabilitation, as well as in the home environment. The design of our standing devices enables independent, stable and safe standing. An electrical belt retractor is used to move the patient from the wheelchair to a secured standing position.

4.2 Indication

People with a height of 140 to 200 cm with various physical and/or mental limitations are enabled to stand safely, firmly and without fatigue, which is necessary every day. Daily standing stabilizes the circulation, prevents bone decalcification (osteoporosis), corrects contractures of the hip and knee joints, stimulates bowel and bladder functions as well as respiration and kidney activity. The targeted promotion of head control is optimized by the precise adjustment options at all levels.

4.3 Contraindication

The following patients are not permitted to use the standing device:

- In case of decubitus especially in the hip and leg area
- In case of extreme deformation and non-weight-bearing capacity of the lower extremities
- In case of massive cardiovascular problems
- In case of severe dizziness, which makes it impossible to stand up independently or partially independently
- In case of severe anxiety

4.4 Equipment features

The Levacare standing device has the following features:

- Table cutout with ergonomic abdominal pad
- Foot plate with heel cups and toe fixation
- Double swivel castors with total locking
- Table height adjustment with gas spring support
- Knee pads adjustable in height and depth, as well as horizontal and angular adjustment
- Stainless steel handrails parallel height adjustable
- Belt retraction system available mechanical and electrical
- Easy maintenance and repair of the electrical or mechanical belt retraction system

4.5 Scope of delivery

- Pre-installed stable base frame with bar rails adjustable from 84 to 107 cm
- Support table with height adjustment from 100 to 130 cm
- Height and depth adjustable knee pads
- Heel rests with velcro strap fixations

Optional:

- Stand-up aid with electric belt automatic

5 Assembly information

5.1 Basic information for assembly

The standing frame should only be assembled by authorised personnel. The fuse on the installation side must not exceed 16A.

Before connecting the device, please ensure that the voltage and frequency of your mains supply correspond to the specifications on the type plate.

Ensure that the standing frame is placed on a level surface when selecting its location. Ensure a suitable floor covering if the standing frame has to be moved frequently. Carpets, rugs and loosely laid floor coverings can be damaged or make it difficult to push the standing frame.

Connect the mains plug firmly to the mains socket. When doing so, lay the mains connection cable on the floor. When doing so, ensure that the standing frame (especially when moving) is not standing with its castors on the cable. The cable must not be fed through the mechanics of the standing frame! (danger of crushing)

Damage to the electrical mains cable by running over it or clamping it can have fatal consequences.

6 Operation

All presettings may only be made on the empty standing frame without the patient standing in it. The settings must then be checked again with the patient standing in the aid and corrected if necessary. An assistant is required for the adjustment work. Due to the large number of setting options for the standing frame, some of the settings described may not be possible for some users or in certain constellations.

After adjustment work, check all connections for tightness!

The padding of the table must never be set below the top edge of the basin!

6.1 Body height

The adjustment of the body height is made via the extension height of the table tubes. The setting is correctly selected when the therapy table is at arm's length height with the patient standing or slightly above this required (cf. Figure 2). The table height adjustment is supported by two built-in gas springs.

Make the clamping lever loosen at the tube and pull it out if necessary (cf. Figure 3)

While using the standing frame the clamping lever must always be tightened.

Figure 2: Adjustment of the table height

Figure 3: Clamping lever for table height adjustment

6.2 Lateral pads

An assistant is necessary for the use of the lateral pads!

The height of the lateral pad is adjusted by loosening the clamping levers on the main tubes and the depth is adjusted by loosening the clamping levers on the rear cross tube. Adjustment to the width is made by means of catches on the rear cross tube. The pads are positioned above the buttocks in the area of the lumbar vertebrae with the patient standing (see Figure 4). They are intended to have a stabilizing effect in conjunction with the heel holder and knee pad.

While using the standing frame the clamping lever must always be tightened!

Figure 4: Positioning the lateral pads

6.3 Foot and heel holder

Adjust the heel holders by repositioning them on the perforated grid plate so that the user has a firm hold in them (4 positions are possible).

The patient can step in by opening the Velcro straps on the belts. The straps are closed tightly over the instep. The Velcro straps serve as toe fixation. (cf. Figure 5)

Check the firm hold of the foot!

Figure 5: Adjustment of heel cups and toe fixation

6.4 Knee pads

The knee pads are adjusted in height or depth by loosening the vertical or horizontal tube snaps. The knee pads are correctly adjusted when they rest against the tibia below the knee. The knee pads should be preset first and corrected if necessary when the patient is standing. They stabilize the patient's stance in conjunction with the heel holder and back pad. (cf. Figure 6)

Figure 6: Adjustment of the knee pads

6.5 Thoracic Pad

6.5.1 Setting the Thoracic Pad

The Thoracic Pad provides support to the centre of the back at shoulder blade height. To fit the pad, turn the knob on the support bracket through one full turn and pull it out. Slide the Thoracic pad tube into the bracket with the holes facing out.

Lower the tube into the bracket and release the knob once the thoracic support is at the desired height. Let go of the knob so that the indexing pin built into the knob locates into the nearest hole in the tube. Tightening the knob clamps the Thoracic pad in position.

6.5.2 Thoracic Lateral Pads

The Thoracic Pad is complimented with thoracic lateral pads. Two width adjustable pads can be fitted onto the thoracic pad to fit on either side of the chest at axilla (armpit) height. The height of the pads is dependent on the height setting of the thoracic pad and should be set so that there is a three finger gap between the top of the thoracic lateral pads and the axilla. See previous instructions for setting thoracic pad height.

6.5.3 Adjusting the Width Setting of the Thoracic Lateral Pads

The thoracic supports can be individually adjusted laterally to centralise the chest (thorax). Their position is held with indexing plungers. Pull the plungers out and slide the thoracic supports along the horizontal tube to the desired position. Release the plungers into the nearest set of holes and ensure they have popped into position. Aim to position the supports symmetrically.

6.5.4 Fitting and adjusting the Head Support

The Head Support provides support to the head at occiput height (base of skull, same height as top of ears). To fit the Head Support, pull out the indexing plunger on the thoracic pad. Slide the Head Support tube into the bracket with the holes facing out. Once at the desired height, release the indexing plunger so that it locates positively into the nearest hole in the headrest tube.

The Head support depth should be adjusted so that the head is supported but not pushed forwards. To adjust the depth of the head support, pull up the plunger on the top of the head support and move

the support in or out to the desired position. Once in the desired position, release the plunger and ensure it registers with one of the holes on the head support.

Figure 7: Thoracic Pad

6.6 Thorax-Support

The thorax-support supports the patient while standing. After centering and aligning the thorax-support on the table, it can be attached to the table of the standing device with two screw clamps.


Figure 8: Thorax-support

6.7 General instructions for use when standing up

Fasten the patient belt before approaching the standing frame.

The patient belt should be fastened so that its lower edge is clamped between the buttocks and the wheelchair seat. Then close the abdominal belt. (cf. Figure 9)

Swing the leg rests or footrest to the side. Place the feet in the foot guides. The heel cups can be removed for this purpose. Move the wheelchair so far that the front wheels touch the recesses provided for this purpose in the footplate. Securing the wheelchair against rolling away. (cf. Figure 5)

Figure 9: Position of the patient support belt in the buttocks when the patient is in the wheelchair

After the electric retractable belt (cf. Figure 11) has been connected to the power supply using the mains plug, the device is ready for operation. The retractor can be operated with the hand control (cf. Figure 10) or optionally with the suck/blow tube. The drive is automatically switched off in the two end bearings of the webbing. This prevents the patient from becoming trapped under any circumstances

Figure 10: Hand control

The additional forced shutdown is triggered as long as the shutdown button (cf. Figure 11) is pressed. The unit can only be operated again after it has been released. Each upward or downward movement of the belt straps is effected by pressing the respective key of the manual switch. The keys are marked with arrows according to their function. If suction is applied to the mouthpiece of the hose, the webbing straps are moved up and down when blowing at the mouthpiece. The air hose must always be laid and stored in such a way that unintentional functions by the user, assistant or third parties are excluded. Malfunctions can occur if the air hose is kinked or crushed. If the air switch is not required, the air hose can be pulled off the connector.


Figure 11: Electric retractable belt

6.8 Operating instructions

Before using the standing frame with electric retractable belt, an individual adjustment is required. In particular, a correct adjustment of the table height and the switch-off button of the forced switch-off is important for the safe function of all devices. The table height must be adjusted by loosening the clamping levers in such a way that the padding of the table top is located between the sternum and the upper edge of the pelvis in the standing position according to personal requirements. (see chapter 6.1)

The padding of the table top must never be set below the top edge of the basin!

After adjustment, tighten the clamping levers again.

To set the switch-off button of the forced switch-off, loosen the lock nut and turn the switch-off button on the thread of the switch rod. The switch-off button is set correctly if the motor is not yet switched off up to the desired standing position, but the switch-off is imminent.

After the setting has been made, the shut-off button must be secured again with the lock nut.

When replacing or changing the belt retractor, make sure that no person is standing by the unit when this is being done. Furthermore, the device must be disconnected from the power supply.

Now remove the two locking screws on the two guide rods. Approach the unit from behind and support the belt retractor with your forearms to prevent it from slipping uncontrollably when you pull the locking bolts. You can now let the belt retractor slide down the retaining rods and remove it.

Mounting is carried out in reverse order. The unit is correctly mounted when the indexing bolts are correctly engaged and the locking screws are attached

While the forced switch-off is actuated, all functions of the belt retractor are switched off. Further operation is only possible again when the switch-off button is released, if necessary press it back by hand. When unwinding the belts, make sure that they are tensioned to prevent the belts from becoming entangled.

6.9 Behavior in case of power failure

If the power fails during use of the electric belt retractor, then an assistant must be called by the user to help the patient back into the wheelchair.

The following description assumes that the assistant is right-handed. If the patient is left-handed, proceed the other way round.

First, the wheelchair is turned around the left front wheel just enough so that the assistant can stand behind the patient, but in the process it is possible to place the patient in the wheelchair. The assistant stands behind the user, opens the Velcro fastener, grasps the patient's chest with the right arm, passing under the patient's right armpit, and opens the left belt buckle of the patient restraint belt with the left hand (if necessary, relieve the patient restraint belt by applying a little pressure over the body buckle). The assistant's left hand now grips under the patient's buttocks to help him slide back into the wheelchair.

6.10 Operating the brake castors

The standing frame should always be braked at the place of installation with the help of the castor brake.

Figure 12: Exemplary castors in braked and unbraked condition

Standing training with non-braked castors can lead to patient injury.

7 Ambient conditions

7.1 Storage conditions

Storage conditions	min. -10 °C max. +50 °C	
Relative humidity	min. 20 % max. 80 %	
Air pressure (at altitude ≤ 3000 m)	min. 700 hPa max. 1060 hPa	

Table 3: Storage conditions

7.2 Operating conditions

Operating conditions	min. +5 °C max. +40 °C	
Relative humidity	min. 20 % max. 80 %	
Air pressure (at altitude ≤ 3000 m)	min. 700 hPa max. 1060 hPa	

Table 4: Operating conditions

8 Technical data

(Subject to change without notice)


Designation	KR-100-B	KR-111-B
Nominal voltage		~230 V/50Hz
Rated power		690 VA
Device type B according to IEC 601-1		
Protection class		
Sound power level		63 dB(A)
IP protection class Hand switch, control box, motor		IP X4
IP protection class standing frame mechanism		IP X1
Duty cycle switch on duration 10%		maximal 6 Minuten/Stunde
Max. patient weight	120 kg	120 kg
Patient size	140 - 200 cm	140 - 200 cm
Weight	ca. 52 kg	ca. 65 kg
Width	77 cm	77 cm
Length	118 cm	118 cm
Table height	100 - 130 cm	100 - 130 cm

Table 5: Technical data

Reparations may only be carried out by ISKO qualified personnel or by persons authorized and trained by ISKO with comprehensive product knowledge. In case of non-compliance with this provision, any warranty and liability claims will be rejected.

9 Used materials

The medical device is manufactured as a welded tubular steel construction. The surfaces are powder coated or galvanized. All wooden parts are either laminated or lacquered. The surfaces of this product are unthinkable for the skin from the point of view of health.

10 Service and care

All household cleaners without ammonia and scouring agents are permissible for cleaning the tube parts, the lying surface and the wooden parts with a damp cloth. Solvents (e.g. nitro) destroy the coating of the tubes!

Mechanical cleaning (e.g. scraping, sanding) or jet cleaning of the standing frame is not permitted. All pivots of the moving parts, including the bearing eyes on the adjustment device, are provided with maintenance-free slide bearings and must not be oiled or greased.

11 Service life of the product

With an expected average degree of use in home care, the service life of the standing frame is 10 years. Lack of maintenance and excessive stress on the product can significantly reduce the service life of the standing frame. The expected service life in professional nursing home operation is 7 years.

12 Disinfection

- In order to ensure that the standing frame functions properly, each ISKO standing frame should be cleaned, disinfected and checked after each use so that it can be used again immediately.
- Improper cleaning/disinfection of the standing frame can cause hazards.
- Depending on the degree of soiling, we recommend cleaning the standing frame with a damp cloth or similar.
- For wipe and spray disinfection, disinfectants in their intended concentration can be used. (see manufacturer's instructions)
- The dilution ratio recommended by the manufacturers in the respective instructions for use must be used.

Solvents are not permitted.

Abrasives or scouring sponges must not be used.

12.1 Specifications of detergents and disinfectants

- The working solutions should normally be used freshly prepared.
- The concentrations given should not be exceeded or fallen below.
- They must not contain corrosive or caustic components.
- They must not contain any substances that alter the surface structure or the adhesion properties of the materials.
- Lubricants must not be attacked by cleaning and disinfecting agents.

Under no circumstances should soap or washing-active substances be added to the disinfectant. In the case of products containing alcohol, there is a risk of explosion and fire when applied over large areas.

The use of unsuitable detergents and disinfectants can cause damage to the surface coating for which ISKO KOCH GmbH cannot be held liable.

13 Operational faults and solutions

	Fault	Measure
1	None of the motors respond to the switch actuation.	Check plug connection between the hand control and control box. Check plug connection between the motor cables and control box. Check the plug connection of the power cord.
2	A motor does not respond to the switch actuation	Check plug connection between the motor cable and control box. Check plug connection between the hand control and control box.

Table 6: Operational faults and solutions

For issues which cannot be rectified using the aforementioned instructions; any changes, new settings or repairs to the standing frame may only be implemented by the manufacturer directly, or by a workshop authorised by the manufacturer.

14 Recommended accessories

Article	Order number
Support bar with hip pads	KR-620-X
Thoracic & Head Support	KR-631-X
Thorax-Support 15cm	KR-744-X
Thorax-Support 25cm	KR-745-X
Sip & blow control	KR-601-SB
Table plate horizontal adjustable	KR-710-X
Foot plate horizontal adjustable	KR-720-X
Table edge	KR-602-X
Electric belt retractor	KR-601-X
Lift Support harness, Size S, Standard	KR-625-X
Lift Support harness, Size M, Standard	KR-626-X
Lift Support harness, Size L, Standard	KR-623-X
Lift Support harness, Size S, with lip	KR-608-X
Lift Support harness, Size M, with lip	KR-607-X
Lift Support harness, Size L, with lip	KR-606-X
Lift Support harness, Size XL, with lip	KR-609-X

Table 7: Recommended accessories

15 Maintenance

15.1 Legal basis

The Medical Device Regulation (EU) 2017/745 (MDR) as well as national laws and regulations require operators of medical devices to ensure a safe operating condition of the medical device during the entire period of use.

15.2 Maintenance intervals

As a requirement of the Medical Device Operator Ordinance §4 (Maintenance), a thorough visual inspection (1), a functional test (2) and a current discharge test (3) must be performed in accordance with DIN EN 62353:2015-10 after the medical device has been in operation for at least two years.

(1) During the visual inspection, particular attention must be paid to the following points:

- Tight fit of all screw connections
- Mobility of the pivot points
- Checking the power supply cable for pinching or shearing points
- Check of the strain relief of the power supply line

(2) During the functional test, special attention shall be paid to the following points:

- Function of all electrically operated movements
- Fully extend and retract all motors on the standing frame until they switch off by themselves. (Limit switches in the motors must switch off with an audible click).
- Functionality of the brakes
- Mobility of the triggers
- Check of the hand control

Functional tests and current leakage tests may only be carried out by ISKO specialist personnel or by persons authorized and trained by ISKO with comprehensive product knowledge.

15.3 Spare parts

All spare parts for this medical device must be ordered from IKSO KOCH GmbH, stating the serial number, order number and article number (these can be found on the type plate attached to the medical device).

To ensure that the functional safety and any warranty claims remain valid, only original ISKO KOCH GmbH parts are to be used for the spare parts.

ISKO KOCH GmbH

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95448 Bayreuth

Tel.: +49(0)921/150845-0 (Monday – Thursday 8:00 – 17:00 pm & Friday 8:00 – 16:00 pm)

Fax: +49(0)921/150845-45

E-Mail: info@isko-koch.de

15.4 Notes on documentation

According to the Medical Device Operator Regulation and Medical Device Regulation (EU) 2017/745 (MDR), there is a documentation obligation for:

- Maintenance
- Incidents / near misses

If extraordinary hazards for the product are foreseeable at the installation site of the standing frame (supply line lies on the floor; children playing; pets; ...), the electrical lines in particular must be constantly checked and suitable measures taken to avert hazards.

16 Reuse

Before each reuse of the standing frame, a thorough visual and functional check of all electrically operated functions as well as a current leakage test according to DIN EN 62353:2015-10 must be performed as described under the item Maintenance intervals. The points on service and care & maintenance mentioned in the operating instructions must always be observed when cleaning the standing frame.

17 Disposal

17.1 Disposal of the device

Disposal of the device and accessories, if any, should be carried out in an environmentally friendly manner and in accordance with the legal regulations. Please adhere to the valid waste separation regulations! If there are any uncertainties in this matter, please contact your local municipality or waste disposal company.

17.2 Disposal of the electrical components

*if electrical components are included in the medical device

According to Directive 2012/19/EU - WEEE2, this medical device is classified as an electrical device. All electrical components are free of unauthorized ingredients classified as harmful according to RoHS II Directive 2011/65/EU. In addition, replaced electrical components must be disposed of in accordance with European directives (see Directive 2012/19/EU - WEEE2).

17.3 Disposal of the packaging

The EU Waste Framework Directive 2008/98/EC is decisive for the handling during the disposal of the packaging. Reusable materials must be fed into a recycling cycle in accordance with national regulations.

18 Declaration of Conformity

As the manufacturer, we declare under our sole responsibility that our standing frames complies with the basic requirements of the EC Directive for

Medical devices 2017/745, Annex II

ISKO KOCH GmbH

Egerländer Str. 28

95448 Bayreuth