

# Exaflip

The space-saving wall couch

With  
or without  
electric height  
adjustment

Electrical  
height adjust-  
ment range:  
400 mm

Space-saving  
and versatile

Individually  
configurable



◀ The version without the electric height function can be installed variably on the wall so that the lying height can be selected as required.

▶ When folded up, space is created for other room requirements.



The **Exaflip** wall couches combine flexibility with functionality. Wall-mounted couches are the perfect solution for use as needed, e.g. in first aid rooms in schools, public authorities or companies. The wall couch from ISKMed is individually configurable and offers a wide range of accessories.

## Product details


- | Lying surface dimensions: 70 – 90 x 140 – 200 cm
- | Lying height ranges:
  - without height adjustment: 66 – 86 cm
  - with electrical height adjustment: 66 – 106 cm
- | Pneumatically assisted folding function
- | Hand control with magnetic holder
- | Lockable height function



**Product enquiry**  
**Exaflip**

## Product images



The product images on this page show the electric version of the wall couch. On image No. 3 you can see the variant with fixed height. The final product changes visually due to optional equipment ("Extra equipment/optional" on next page), such as side rails. If you have any questions regarding an order, please do not hesitate to **contact us** .

## Technical data

- | Width: 70 – 90 cm
- | Length: 140 – 200 cm
- | Lying height ranges:
  - 66 – 86 cm (fixed height)
  - 66 – 106 cm (electric couch)
- | Safe working load: 135 kg

## Extra equipment / optional

- | Fixed wall couch with gas spring support foldable upwards (WL-002-0)
- | Electrically operated wall couch foldable upwards (WL-003-0)
- | Side rails (WL-202-0)
- | Battery operation
- | Mounting rail (WL-012-0)

