

# Exasan Duo

The practice couch with head section

Head-down  
position  
possible

Head section  
ratchet  
adjustment:  
-30° - +45°

4 adjustable  
screw feet



◀ The integrated ratchet adjustment offers a wide range of adjustment options for the head section.

▶ The couch can optionally be fitted with a paper roll holder.



The **Exasan Duo** practice couch has an adjustable head section, which enables more specialised examination variants in the treatment room. The stable design and adjustable screw feet provide safety for patients and medical professionals.

## Product details

- | Lying surface: 70 x 200 cm
- | Head section dimensions: 70 x 45 cm
- | Standard height: 50 cm
- | Available with nose slit incl. cover

## Product images



The product images on this page show the standard version of the practice couch. The final product changes visually due to optional equipment ("Extra equipment/optional" on next page), such as a nose slit incl. cover. If you have any questions regarding an order, please do not hesitate to [contact us](#) .

## Extra equipment / optional

| Nose slit incl. cover (BL-010-N)

| Paper roll holder (BL-151-0)



**Product enquiry**  
**Exasan Duo**

**ISKO Koch GmbH** · Egerländer Str. 28 · 95448 Bayreuth (GER)  
Tel +49 921 150845-0 · [info@iskomed.de](mailto:info@iskomed.de) · [www.iskomed.de](http://www.iskomed.de)