

Levacare / Levacare Pro

The standing device for rehab and your home

Height
adjustable
belt system

Pneumatically
assisted
table height
adjustment

Stainless steel
handrails
height-adjust.
in parallel

Max.
patient
weight
120 kg



◀ The multifunctional knee pads are included as standard with all standing devices.

They can be adjusted horizontally, vertically and in depth and have an integrated angle adjustment that adapts to the leg/ knee position during use.





Standing devices are among the most important devices in rehabilitation and care. The stable design of the **Levacare** gives patients a firm, secure hold during their daily standing training. The standard version of the Levacare is supplied without a belt system, while the Pro version has an electric stand-up device.

Product details

- | Table cut-out with ergonomic belly pad
- | Foot plate with heel cups and toe fixation
- | Lift support harness (size S, M or L)
- | Sturdy tubular steel construction in blue or brown
- | Wooden elements made of multiplex with high-strength HPL coating
- | 4 Double swivel castors with total locking
- | Integrated suction-blow control on the Levacare Pro

Product images



The product images on this page show the Levacare Pro variant of the standing device. In the Levacare variant, the electric belt system is omitted. The final product also changes visually due to optional equipment ("Extra equipment/optional" on next page), such as a 4-part table edge made of beech. If you have any questions regarding an order, please do not hesitate to contact us .

Technical data

- | Total length: 118 cm
- | Table width: 66 cm
- | Total width: 77 cm
- | Table length: 59 cm
- | Tie bar height: 84 – 107 cm
- | Patient weight: 120 kg
- | Tie bar width: 48 cm
- | Min./Max. patient size: 120 – 220 cm
(smaller than 140 cm, a height-adjustable
foot plate KR-724-X must be ordered)
- | Table height: 100 – 130 cm

Extra equipment / optional

- | Levacare (KR-100-B) without belt system
- | Levacare Pro (KR-111-B) with electric belt system
- | Lift Support harness „S“ (KR-608-X)
- | Lift Support harness „M“ (KR-607-X)
- | Lift Support harness „L“ (KR-606-X)
- | Lift Support harness with rubber (KR-604-NO)
- | Lift Support harness with sitting assist
- | Lift Support harness with bigger bottom
- | Support bar with hip pads (KR-620-X)
- | Upper body support incl. headrest (KR-631-X)
- | Table top horizontally movable (KR-710-X)
- | Foot plate horizontally movable (KR-720-X)
- | Foot plate horizontally adjustable & height adjustable (KR-724-X)
- | Breast pad 15 cm (KR-744-X)
- | Breast pad 25 cm (KR-745-X)
- | Knee pads extra soft (KR-801-X)
- | Heel Velcro (KR-805-X)
- | Table edge beech, 4-piece (KR-602-X)
- | Foot plate elevation