

Levaflex

The ergonomic standing device

standing devices



Swivelling rear wheels

Synchronised operation of belt and tilt motor

Pneumatically assisted table height adjustment

Max. patient weight 120 kg



< The tilt function makes it easier to get up and put down and is less stressful for the patient.

> The tilt and automatic seat belt functions can be switched synchronously or separately.




The **Levaflex** is the standing device with the highest functionality. The arm rails, knee support and footplate are tilted towards the patient before standing up or sitting down, making the process easier. The patient's knees and behind are under less strain thanks to the tilt function, with the usual stability and safety.

Product details

- | Ergonomically adjustable knee guide
- | Table cut-out with ergonomic belly pad
- | Foot plate with heel cups and toe fixation
- | Electric belt system
- | Lift support harness (size S, M or L)
- | Sturdy tubular steel construction in blue or silver
- | Wooden elements made of multiplex with high-strength HPL coating
- | Swivelling rear wheels
- | 4 double swivel castors with total locking

Product images



The product images on this page show the standard version of the standing device. The final product changes optically due to optional equipment ("Extra equipment/optional" on next page), such as a 4-piece table edge made of beech. If you have any questions regarding an order, please do not hesitate to contact us .

Technical data

- | Total length: 100 – 103 cm
- | Table width: 64 cm
- | Total width: 73 – 78 cm
- | Table length: 58 cm
- | Tie bar height: 80 – 100 cm
- | Patient weight: 120 kg
- | Tie bar width: 49 cm
- | Min./Max. patient size: 140 – 220 cm
- | Table height: 100 – 132 cm

Extra equipment / optional

- | Levaflex (KR-214-B) frame colour blue
- | Levaflex (KR-214-S) frame colour silver
- | Lift Support harness „S“ (KR-608-X)
- | Lift Support harness „M“ (KR-607-X)
- | Lift Support harness „L“ (KR-606-X)
- | Lift Support harness with rubber (KR-604-NO)
- | Lift Support harness with sitting assist
- | Lift Support harness with bigger bottom
- | Support bar with hip pads (KR-620-X)
- | Upper body support incl. headrest (KR-631-X)
- | Table top horizontally movable (KR-710-X)
- | Breast pad 15 cm (KR-744-X)
- | Breast pad 25 cm (KR-745-X)
- | Knee pads extra soft (KR-801-X)
- | Heel Velcro (KR-805-X)
- | Battery box (KR-703-X)
- | Table edge beech, 4-piece (KR-602-X)
- | Foot plate elevation