

# Levatrain

The functional and inexpensive standing device

standing devices



Electric belt system with premium linear drive

Pneumatically assisted table height adjustment

Stainless steel handrails height-adjust. in parallel

Max. patient weight 120 kg



- < The integrated table height adjustment is very smooth-running thanks to gas spring support.
- > The standard swivelling rear wheels make it easier to drive the wheelchair into the standing device. This is particularly advantageous for wheelchairs with wider seats.



The **Levatrain** standing device has an integrated electric belt system and makes it easier to stand up while providing a stable support. The trapezoidal base frame and the swivelling rear wheels also allow the standing device to be used by patients with wider wheelchairs.


## Product details

- | Table cut-out with ergonomic belly pad
- | Foot plate with heel cups and toe fixation
- | Electric belt system
- | Lift support harness (size S, M or L)
- | Sturdy tubular steel construction in blue or silver
- | Widenable access area
- | Wooden elements made of multiplex with high-strength HPL coating
- | Swivelling rear wheels
- | 4 double swivel castors with total locking



## Product images



The product images on this page show the standard version of the standing device. The final product changes optically due to optional equipment ("Extra equipment/optional" on next page), such as a 4-piece table edge made of beech. If you have any questions regarding an order, please do not hesitate to [contact us](#) .

## Technical data

- | Total length: 100 – 103 cm
- | Table width: 64 cm
- | Total width: 73 – 78 cm
- | Table length: 58 cm
- | Tie bar height: 80 – 100 cm
- | Patient weight: 120 kg
- | Tie bar width: 49 cm
- | Min./Max. patient size: 140 – 220 cm
- | Table height: 100 – 132 cm

## Extra equipment / optional

- | Lift Support harness „S“ (KR-608-X)
- | Lift Support harness „M“ (KR-607-X)
- | Lift Support harness „L“ (KR-606-X)
- | Lift Support harness with rubber (KR-604-NO)
- | Lift Support harness with sitting assist
- | Lift Support harness with bigger bottom
- | Support bar with hip pads (KR-620-X)
- | Upper body support incl. headrest (KR-631-X)
- | Table top horizontally movable (KR-710-X)
- | Breast pad 15 cm (KR-744-X)
- | Breast pad 25 cm (KR-745-X)
- | Knee pads extra soft (KR-801-X)
- | Heel Velcro (KR-805-X)
- | Battery box (KR-703-X)
- | Table edge beech, 4-piece (KR-602-X)
- | Foot plate elevation