

# Multidorm Mobil

The nursing bed with electric folding function

Max.  
patient  
weight  
135 kg

Safe  
working load  
175 kg

Dimensions:  
90 x 200 cm

Height  
adjustment  
range  
40 – 80 cm






◀ After the automatic unfolding process, the Multidorm Mobil is a fully-fledged nursing bed that complies with the requirements of the ISO 60601-2-52 standard.


The **Multidorm Mobil** is perfect for quick use in home care. It can be folded and unfolded electrically at the touch of a button and can be stored in a space-saving way. Thanks to the patented side rail solution, no tools are required to set up and dismantle the bed.

## Product details

- | Mobile single bed with 90 cm width
- | Trendelenburg and anti-Trendelenburg positioning possible
- | Individually braked and very robust castor system
- | Particularly fast assembly and disassembly due to integrated folding motor
- | Integrated holding brackets for transporting the side rails
- | Patented side rail pick-up system
- | Easy to transport and store thanks to small dimensions
- | Can be unfolded by one person alone

## Product images



The product images on this page show the standard version of the care bed. The final product changes visually due to optional equipment ("Extra equipment/optional" on next page), such as a side rail cover made of imitation leather or fabric. If you have any questions regarding an order, please do not hesitate to [contact us](#) .

## Technical data

- | Mattress width: 90 cm
- | Mattress length: 200 cm
- | Bed height: 40 – 80 cm
- | Max. patient weight: 135 kg
- | Safe working load: 175 kg

## Extra equipment / optional

- | Side rail extension (NS-058-0)
- | Side rail cover imitation leather (NS-049-1)
- | Side rail cover fabric (NS-049-0)