

Rotadorm Max / Max Low

The stand-up bed for obese patients


Max.
patient
weight
180 kg

Safe
working load
220 kg

3 lengths
adjustable:
200, 210 and
220 cm

Height
adjustment
range
34 – 78 cm or
40 – 80 cm

Direction
of rotation
possible on
both sides



◀ The independently controllable neck support ensures greater sitting and lying comfort and can provide ergonomic support when getting up (left image).

◀ Beech-look side panelling creates a cosy, homely ambience.

The Maxi version of our **Rotadorm bed series** offers obese patients stability in everyday life, whether sleeping or getting up. The Low version is ideal for people under 165 cm tall.






Product details

- | 100 cm wide single bed, adjustable to three lengths
- | Height, head section, neck section, thigh support, lower leg support and turning mechanism individually electronically adjustable
- | 2-section heavy-duty mattress with incontinence cover
- | Patient lifting pole with handle
- | Removable patient handles in frame colour
- | Programming and patient hand control
- | Optionally with battery system for self-sufficient use of all functions
- | Stable scissor lift system for height adjustment
- | Individually braked double swivel castors
- | Emergency release for rotating mechanism for Rotadorm Max Low




Product enquiry
Rotadorm Max / Max Low

Product images



Max Low



Max Low



Max Low



Max Low



Max Low

The product images on this page show both the version of the Rotadorm Max with a height adjustment range of 40 – 80 cm, and the version Max Low (images are marked accordingly) with a height adjustment range of 34 – 78 cm. In addition, in the images you can see the version of Rotadorm Max with mattress both with and without white cover (the cover is optionally available). The final product changes visually due to optional equipment ("Extra equipment/optional" on next page), such as the battery pack. If you have any questions regarding an order, please do not hesitate to [contact us](#) ✉.

Technical data

- | Mattress width: 100 cm
- | Mattress length: 200 / 210 / 220 cm
- | Bed height: 40 – 80 / 34 – 78 cm
- | Max. patient weight: 180 kg
- | Safe working load: 220 kg

Extra equipment / optional

- | Rotadorm Max (bed height 40 – 80 cm)
- | Rotadorm Max Low (bed height 34 – 78 cm)
- | Lying surface dimension: 100 x 200 cm (SI-021-1)
- | Lying surface dimension: 100 x 210 cm (SI-021-21)
- | Lying surface dimension: 100 x 220 cm (SI-021-22)
- | Turningside: left
- | Turningside: right
- | **Cover for handle receptacle (SI-011.80.725)**
enables getting out of bed to the side
- | **Battery pack with charging circuit (SI-011.85.050)**
safety equipment for power failure – all electrical functions can thus be returned to the initial position even without mains voltage

