

Rotadorm Vario

The installation frame with stand-up function

Max.
patient
weight
up to 135 kg

Safe
working load
up to 175 kg

Height
adjustment
range
40 - 80 cm

For bed frames
with a width
of 90 cm

Direction
of rotation
possible on
both sides



◀ The stand-up function of the Vario version of our Rotadorm bed series can also be customised to the patient and saved.

▶ With a special installation frame placed in any bed, a care bed can be created that perfectly matches the interior design style.




The **Rotadorm Vario** is a special installation frame which, in addition to the usual comfortable functions of ISKMed nursing beds, also has a turn/sit and stand-up function. The frame fits perfectly into an existing bed frame and offers patients and carers every comfort in everyday care. The Rotadorm Vario can be used in a single or double bed.

Product details

- | High-quality 6-motor nursing bed system
- | Length individually adjustable: 200 and 220 cm
- | 2-section special mattress with integrated brackets
- | 5-section lying surface with springwood slats
- | Variable adjustable mattress support brackets
- | Patient lifting pole with ergonomic handle
- | Programming and patient hand control

Product images



The product images on this page show the optional version with castors. The standard version of the special installation frame is supplied with base columns. The final product changes visually due to optional equipment ("Extra equipment/optional" on next page), such as a pull handle. If you have any questions regarding an order, please do not hesitate to **contact us** .

Technical data

- | Mattress width: 85 cm
- | Mattress length: 200 – 220 cm
- | Bed height: 40 – 80 cm
- | Max. patient weight: 135 kg
- | Safe working load: 175 kg

Extra equipment / optional

- | Turningside: left
- | Turningside: right
- | **Pillow clips (SI-038-0)**
holding the pillow in position when the headboard is up,
to make it possible to stand up comfortably
- | **Pull handle (SI-011.85.307)**
serves as a positioning aid in the bed, is mounted on the headboard
- | **Cover for handle receptacle (SI-011.80.725)**
enables getting out of bed from the side
- | **Table (SI-132-0)**
serves as storage option during sitting position, mounting on handhold
over wire frame, dimensions 365 x 265 x 10 mm
- | **Transport rack (TRANSP-08)**
to assist with transport, deposit € 300
- | **Battery pack with charging circuit (SI-011.85.050)**
safety equipment for power failure – all electrical functions can thus be returned
to the initial position even without mains voltage