

Belladorm Duo

The partner nursing bed

Max.
patient
weight
2 x 135 kg

Maximum
total load
350 kg

Dimensions:
180 x 200 cm

Height
adjustment
range
36 - 76 cm



◀ With a variety of adjustment functions for both halves of the bed, both partners enjoy optimum comfort. The leg elevation designed with a step-by-step locking mechanism is also possible at the touch of a button (left image).

◀ Two separate handsets ensure individual freedom of movement.

Belladorm Duo is the care bed for couples who want to lie and sleep together in a stable and comfortable way. Among other things, the care bed is equipped with two handsets: one controls the head and foot area on one side of the bed. The other handset adjusts the height as well as the head and foot area and also has an auto-contour function for the second side of the bed.


Product details

- | Double bed with individual lying surfaces
- | 2 x 4-section lying surface with stable metal grid support
- | Patient lifting pole with handle
- | 4-way manual switch with single function lock
- | 2-way manual switch with single function lock for the partner bed
- | Stabilisation system for one-sided loading
- | High-quality Linak motor system



Product images



The product images on this page show the standard version of the nursing bed. The final product changes visually due to optional equipment ("Extra equipment/optional" on next page), such as side rail coverings made of imitation leather or fabric. If you have any questions regarding an order, please do not hesitate to [contact us](#) .

Technical Data

- | Mattress width: 180 cm
- | Mattress length: 200 cm
- | Bed height: 36 – 76 cm
- | Max. patient weight: 2 x 135 kg
- | Safe working load: 350 kg

Extra equipment / optional

- | Easy EX stand-up aid (NS-035-0)
- | Easy EX light stand-up aid (NS-035-LC)
- | Side rail cover imitation leather (NS-049-1)
- | Side rail cover fabric (NS-049-0)
- | Patient lifting pole 3-legged (NS-261-0)