

# Multidorm Max

The most stable nursing and stand-up bed

Max.  
patient  
weight  
260 kg

Safe  
working load  
325 kg

Dimensions:  
120 x 220 cm

Height  
adjustment  
range  
43 - 80 cm



special beds



◀ Thanks to its multifunctional mobility, the Multidorm Max is also ideal for use as a stand-up bed in care settings.

The high stability of the Multidorm Max is guaranteed in all adjustable positions.

The **Multidorm Max** has an extremely stable design with a lifting lever system and a larger lying surface, making it particularly suitable for plus-size patients. It offers the best possible support for patients and carers in the heavy-duty care sector. After removing the foot section, it is also possible to stand up in a forward position over the end of the bed.

## Product details

- | Single bed with extra large lying surface
- | Stable patient double lifting pole with adjustable grab handle
- | Stainless steel impact protection made to order
- | Integrated support feet for a safe foot-lowering or standing-up process
- | Trendelenburg and anti-Trendelenburg positioning possible
- | Central braking system with high-quality hospital castors



## Product images



## Technical data

- | Mattress width: 120 cm
- | Mattress length: 220 cm
- | Bed height: 43 – 80 cm
- | Max. patient weight: 260 kg
- | Safe working load: 325 kg

## Extra equipment / optional

Infrared control