

Exadoc Endo

Examination couch for endoscopy

Ergonomically adjustable headrest

Lowerable stainless steel side rails

Central brake system

Stainless steel impact protection / edge protection



◀ The headrest was specially designed for patients lying on their side. It ensures an ergonomic head position during treatment.

▶ Two premium Li-ion battery modules paired with the wall charging system ensure smooth practical operation.



The **Exadoc Endo** was developed specifically for endoscopic examinations and treatments in collaboration with recognized medical experts in the field. The focus was on the functionality of the head section and the fall protection. The couch offers maximum ergonomics and safety for the patient as well as for the attending physician and his staff.


Product details

- | Lying surface: 70 x 200 cm
- | Under-couch clearance: 14,5 cm (area of lifting columns: 10 cm)
- | Head section: 70 x 40 cm
- | Height adjustment: 55 - 95 cm
- | Safe working load: 210 kg
- | Operation with hand switch on the lying surface
- | Integrated safety through automatic locking of electrical functions
- | Unique patient safety system
- | Cordless due to Li-Ion battery system
- | Frame color: RAL 7035 light grey
- | Imitation leather cover freely selectable
- | 100 mm clinic double castors with central brake



Product images



The product images on this page show the standard version of the examination couch. If you have any questions about available colors and properties of the leatherette cover, please do not hesitate to [contact us](#) .

Note on colour representation: The product images shown on this page are for illustrative purposes only. Colour variations between the screen display and the actual product are possible due to different display technologies and individual device settings. This applies in particular to the metal colour RAL 7035 shown. Only the original RAL colour tone according to the manufacturer's specifications is authoritative.

Technical data

- | Width: 70 cm
- | Length: 200 cm
- | Height: 55 - 95 cm
- | Safe working load: 210 kg

Extra equipment / optional

--