

Exadoc Cardio

The ECG couch

**Cardio
cut-out incl.**
padded insert

Head section
electrically
adjustable:
0° - +35°

**Head section
dimensions:**
80 x 70 cm
Cut-out:
30 x 30 cm

**Central
brake system**



◀ The cut-out in the pad makes it easier to perform echocardiography in the lateral position.

▶ The robust central brake system ensures that the table can be positioned quickly and effectively during operation.



The **Exadoc Cardio** examination couch has been specially designed to meet the needs of cardiological practices. The head section has been extended and provided with an echo cut-out on the left side at heart level. Thanks to the large height adjustment range, the medical team can perform the examinations in an ergonomically adapted manner.

Product details

- | Lying surface: 80 x 200 cm
- | Under-couch clearance: 13 cm
- | Height adjustment: 44 – 95 cm with integrated safety lock
- | Lifting and propulsion system based on the principle of integrated safety
- | Safe working load: 210 kg
- | Individually lockable electrical functions
- | Stable chassis with Linak actuators
- | 100 mm clinic double roller system with central brake
- | Magnetic hand control



Product images



The currently delivered model of the ECG (electrocardiogram) couch has slightly modified elements, such as the hand switch and the motor system (control box and motor). Thus, the new version differs visually from the product images shown here. The final product changes visually due to optional equipment ("Extra equipment/optional" on next page), such as side rails. If you have any questions regarding an order, please do not hesitate to **contact us** .

Note on colour representation: The product images shown on this page are for illustrative purposes only. Colour variations between the screen display and the actual product are possible due to different display technologies and individual device settings. This applies in particular to the metal colour RAL 7035 shown. Only the original RAL colour tone according to the manufacturer's specifications is authoritative.

Technical data

- | Width: 80 cm
- | Length: 200 cm
- | Height: 44 – 95 cm
- | Safe working load: 210 kg

Extra equipment / optional

- | Additional foot switch for adjusting the working height (BL-101-1)
- | Side rail for treatment couch, chromed (1 pair) (BL-201-0)
- | Leg section with ratchet adjustment (BL-999-3)
Foot section mechanically adjustable up to +30° with ratchet adjuster
- | Paper roll holder (BL-151-0)
- | Roller lifting system (BL-501-0)

