

Levatrain

The functional and inexpensive standing frame



Electric belt system with premium linear drive

Pneumatically assisted table height adjustment

Stainless steel handrails height-adjust. in parallel

Max. patient weight 120 kg

standing frames



◀ The integrated table height adjustment is very smooth-running thanks to gas spring support.

▶ The standard pivoting rear wheels make it easier to drive the wheelchair into the standing frame. This is particularly advantageous for wheelchairs with wider seats.



The **Levatrain** standing frame has an integrated electric belt system and makes it easier to stand up while providing a stable support. The trapezoidal base frame and the pivoting rear wheels also allow the standing frame to be used by patients with wider wheelchairs.

Product details

- | Table cut-out with ergonomic abdomen pad
- | Foot plate with heel cups and toe fixation
- | Electric belt system
- | Lift support harness (size S, M or L)
- | Sturdy tubular steel construction in blue
- | Widenable access area
- | Wooden elements made of multiplex with high-strength HPL coating
- | Pivoting rear wheels
- | 4 double swivel castors with total locking



Product images



The product images on this page show the standard version of the standing frame. The final product changes optically due to optional equipment ("Extra equipment/optional" on next page), such as a 4-piece table edge made of beech. If you have any questions regarding an order, please do not hesitate to **contact us** .


Technical data

- | Total length: 100 – 103 cm
- | Table width: 64 cm
- | Total width: 73 – 78 cm
- | Table length: 58 cm
- | Tie bar height: 80 – 100 cm
- | Max. patient weight: 120 kg
- | Tie bar width: 49 cm
- | Min./Max. patient height: 140 – 200 cm
- | Table height: 100 – 132 cm


Extra equipment / optional

- | Lift support harness with seat support (5 cm), size “S” (KR-608-X)
- | Lift support harness with seat support (5 cm), size “M” (KR-607-X)
- | Lift support harness with seat support (5 cm), size “L” (KR-606-X)
- | Seat support of lift support harness: 10 cm (additional charge)
- | Non-slip surface with nubs for lift support harness (additional charge)
- | Sitting assistance for lift support harness (additional charge)
- | Support bar with hip supports (KR-620-X)
- | Upper body support incl. headrest (KR-631-X)
(requires support bar for attachment)
- | Table top horizontally movable (KR-710-X)
- | Chest pad 15 cm (KR-744-X)
- | Chest pad 25 cm (KR-745-X)
- | Knee pads extra soft (KR-801-X)
- | Heel Velcro (KR-805-X)
- | Battery box (KR-703-X)
- | Table edge beech, 4-piece (KR-602-X)
- | Foot plate elevation
- | Special colours available (RAL colour system)


Extra equipment / optional




**Lift support harness
with seat support (5 cm)**




**Lift support harness
with seat support (10 cm)**




**Lift support harness
with rubber**



**Lift support harness
with sitting assist**



**Support bar
with hip supports**




**Upper body support
incl. headrest**



Chest pad 15 cm




Chest pad 25 cm




**Knee pads
extra soft**




Heel Velcro




Foot plate elevation




**Table top
horizontally movable**



**Table edge beech
4-piece**



Battery box



**Special colours available
(RAL colour system)**



**Product enquiry
Levatrain**

ISKO Koch GmbH · Egerländer Str. 28 · 95448 Bayreuth
Tel +49 921 150845-0 · info@iskomed.de · www.iskomed.com