

Exaphysio

The configurable physiotherapy couch

Strictly according to the recommendation of the BfArM

Height adjustment range:
44 – 95 cm

Electric or hydraulic height adjustment

Individually configurable



treatment couches



◀ If required, the head section can not only be raised to +30°, but also lowered to -30°.

▶ With the optional roller lifting system, the couch can be brought to an even more stable state.



The **Exaphysio** treatment couch is the functional all-rounder in the physiotherapy practice. The standard version of the Exaphysio can be ordered with a one-piece lying surface. Thanks to numerous variants in terms of size, design and accessories, therapists can configure the couch to suit their requirements and wishes. The high-quality lifting system offers maximum stability for everyday treatment.



Product details

- | Lying surface dimensions freely selectable: 60 – 90 x 140 – 220 cm
- | Under-couch clearance: 13 cm (height increase possible, please inquire)
- | Height adjustment possible electrically or hydraulically
- | Stainless steel foot pedals on both sides (hydraulic option)
- | Safe working load: 210 kg
- | Lifting and propulsion system based on the principle of integrated safety
- | Individually braked castors



Product enquiry
Exaphysio

Product images



Note on colour representation: The product images shown on this page are for illustrative purposes only. Colour variations between the screen display and the actual product are possible due to different display technologies and individual device settings. This applies in particular to the metal colour RAL 7035 shown. Only the original RAL colour tone according to the manufacturer's specifications is authoritative.

The product images on this page show the version of the physiotherapy couch with electric height adjustment. In addition, the currently delivered model of this couch has slightly modified elements, such as the hand switch and the motor system (control box and motor). Additionally shown on the images are optional features, such as the adjustable head section, the nose slit incl. cover as well as the roller lifting system. The final product will therefore change visually due to additional optional equipment ("Extra equipment/optional" on next page), such as chromed side rails or similar. If you have any questions regarding an order, please do not hesitate to **contact us** .

Technical data

- | Width: 60 – 90 cm
- | Height: 44 – 95 cm
- | Length: 140 – 220 cm
- | Safe working load: 210 kg

Extra equipment / optional

- | Electric height adjustment
- | Hydraulic height adjustment
- | Additional foot switch for adjusting the working height (BL-101-2)
- | Side rail for treatment couch, chromed (1 pair) (BL-201-0)
- | Side rail for treatment couch, chromed (1 piece) (BL-202-0)
- | Nose slit in the head section (SB-187-0)
- | Cover for nose slit (SB-187-1)
- | Head section adjustable with ratchet assembly; +/- 30° (BL-310-0)
- | Head section adjustable with gas spring; +/- 30° (BL-311-0)
(surcharge for ratchet adjustment)
- | Head section electrically adjustable; +/- 30° (surcharge for ratchet adjustment) (BL-312-0)
- | Leg section with ratchet adjustment (BL-999-3)
Foot section mechanically adjustable up to +30° with ratchet adjuster
- | Paper roll holder (BL-151-0)
- | Roller lifting system (BL-501-0)
- | Central brake system
- | Padded/foldable surround - for up to 4 sides

- | Pacificblue (No. 432x1385)
Colour imitation leather cover



- | Mediterranean (No. 208x1752)
Colour imitation leather cover



- | Saffron (No. 208x2388)
Colour imitation leather cover



- | Racingred (No. 207x1173)
Colour imitation leather cover



- | Mercedes Date (No. 207x3115)
Colour imitation leather cover



- | Special colour
Wide range of additional colours