

Belladorm Max

The standard nursing bed for plus-size patients

Max.
patient
weight
280 kg

Safe
working load
up to 320 kg

Dimensions:
120 x 200 cm

Height
adjustment
range
22 – 62 cm or
40 – 80 cm



◀ The sturdy lying surface of the Belladorm Max is specially designed for the high loads in the care of plus-size patients. A specially developed mechanism in the head and foot sections ensures stability even when the patient is sitting on one side of the bed.

▶ The five-way manual switch also offers a Trendelenburg and a comfort function.



The **Belladorm Max** is the perfect solution for patients weighing up to 280 kg. The nursing bed offers variable floor clearance for the use of patient lifts. The Belladorm Max is registered as an aid with health and care insurance companies.

Product details

- | 4-section lying surface with metal grid support
- | Patient lifting pole with grab handle (optional)
- | Stabilization system for lateral loads
- | Electric adjustment for height, head and leg section
- | Hand control with 4 or 5 functions and single function lock
- | Trendelenburg and anti-Trendelenburg positioning possible up to 220 kg patient weight
- | High quality Linak motor system
- | Individually braked castors

Product images



The product images on this page show the version of the standard nursing bed with a height adjustment range of 22 – 62 cm. The final product changes visually due to optional equipment ("Extra equipment/optional" on next page), such as side rail coverings made of imitation leather or fabric. If you have any questions regarding an order, please do not hesitate to **contact** us .

Technical Data

- | Mattress width: 120 cm
- | Mattress length: 200 – 220 cm
- | Height adjustment range: 22 – 62 cm (Max. patient weight: 220 kg / Safe working load: 260 kg)
- | Height adjustment range: 22 – 62 cm (Max. patient weight: 280 kg / Safe working load: 320 kg)
- | Height adjustment range: 40 – 80 cm (Max. patient weight: 220 kg / Safe working load: 260 kg)

Extra equipment / optional

- | **Height adjustment range: 22 – 62 cm**
Max. patient weight: 220 kg
Safe working load: 260 kg
- | **Height adjustment range: 22 – 62 cm**
Max. patient weight: 280 kg
Safe working load: 320 kg
- | **Height adjustment range: 40 – 80 cm**
Max. patient weight: 220 kg
Safe working load: 260 kg
- | **Easy EX stand-up aid (NS-035-O)**
- | **Easy EX light stand-up aid (NS-035-LC)**
- | **Side rail cover imitation leather (NS-049-1)**
- | **Side rail cover fabric (NS-049-O)**
- | **Side rail height extension**
- | **Infrared control**
- | **Patient lifting pole 3-legged with grab handle (NS-261-O)**
- | **Patient lifting pole with grab handle**
provides secure and stable support to make it
easier to sit up or change position in bed



Product enquiry
Belladorm Max

ISKO Koch GmbH · Egerländer Str. 28 · 95448 Bayreuth
Tel +49 921 150845-0 · info@iskomed.de · www.iskomed.com